

Science & Motion Sports

SAM BalanceLab 1.2 for Windows

User Manual

Version 1.2



Manufacturer and Distribution

Science&Motion Sports GmbH
Schaefergasse 4
65428 Ruesselsheim
Germany
Phone +49 (0) 6131 906 4440

info@scienceandmotion.com
www.scienceandmotion.com

Preface

Welcome to the User Manual for the BalanceLab software for analyzing the static and dynamical pressure distribution, using the PDM-S measuring platform.

This User Manual provides a basic understanding for operating the software BalanceLab for analyzing the pressure distribution of static and dynamic forces. It will explain the basic set-up and operation of the PDM-S system, and provide tips on preparation for measuring and data acquisition.

The company Science&Motions Sports GmbH i.G. does not assume any liability whatsoever for injury to personnel or players or damage to the device caused by improper use of the BalanceLab software.

All particulars about the Measurement System in this book were prepared, compiled, and edited with the greatest of care. Nevertheless, we cannot guarantee that mistakes do not appear. We advise the user that Science&Motions Sports GmbH i.G. neither provides a guarantee against nor accepts any legal responsibility or any liability for consequences resulting from incorrect statements.

If you notice any mistakes in this manual, or should you have any suggestions for the improvement, we would be very grateful for giving us a short message any time.

In the interests of continuous product development, the manufacturer reserves the right to carry out improvements to this manual and the product described there in at any time and without any further obligation.

Registered trademarks

Different brand names are referred to in this manual. All these product names are used only for clarity's sake or for editorial reasons and are trade marks belonging to the respective companies.

By using the brand names, the trade marks themselves and the rights of the respective proprietors remain protected.

Science&Motions Sports is a registered trade mark and PDMS denotes a feature of the company, Science&Motions Sports GmbH i.G..

The hardware, software, and company names referred to in this handbook are in general brand names, trademarked, or patented.

Copyrights

This document and extracts taken from it may on no account be duplicated without the explicit consent of Science&Motions Sports GmbH i.G.. The content of this document may on no account be used for purposes that have not been approved. An infringement of the copyright is a punishable offence.

© Copyright Science&Motions Sports GmbH i.G.

All rights reserved.


Conventions

In this user manual the following conventions will be used:

Important warnings are shown as following:

	<p>"WARNING" symbols indicate a potential hazard to the health and safety of the users and/or players. The warnings describe the danger involved and how this can be avoided.</p>
<p>CAUTION</p>	

Notes are shown in this way:

	<p>"NOTE" symbols denote a potential danger that can cause damage to or destroy the device. These NOTE symbols describe the danger involved and how this can be avoided.</p>
<p>NOTE</p>	

This manual is to be kept within easy reach so that the information it contains is available to the user at any time.

Important Security Note:

The surface of the system is designed for indoor usage with sports shoes.



- never use the system if damp or humidity is on the surface or on a player's shoes – there is a danger of slipping or falling for the player in this case !
- never step on the plate with golf shoes equipped with hard spikes – the surface of the plate might get damaged over time!

Science&Motion is currently working on a solution with a covering mat for usage with spikes or outdoor usage. Please contact us for availability date.

Table of contents:


1	INTRODUCTION	6
1.1	Layout of the program.....	6
2	PREPARING THE DATABASE (STEP 1).....	7
2.1	Instructions	7
2.2	General information	7
2.3	Create new players	8
3	ANALYZING THE STATICAL PRESSURE DISTRIBUTION	9
3.1	Instruction for the stance analysis.....	9
3.2	Execution of measurement (step 2)	11
3.3	Preparation for the report in the „View“ (step 3).....	14
3.3.1	Visualization of the measurement data.....	15
3.3.2	Preparation	16
3.4	The report (Step 3).....	17
3.4.1	Explanation of the report.....	17
	The histogram represents the average percentage pressure distribution of the left and right fore and back foot and the percentage pressure distribution between left and right foot.	19
3.4.2	Comparison of two records.....	19
3.4.3	Data Rating.....	20
4	ANALYZING THE DYNAMIC PRESSURE DISTRIBUTION.....	FEHLER! TEXTMARKE NICHT DEFINIERT.
4.1	Instructions for the dynamic analysis	Fehler! Textmarke nicht definiert.
4.2	Execution of measurement (Step 2).....	Fehler! Textmarke nicht definiert.
4.3	Preparation for the report in the “View“ (Step 3)	Fehler! Textmarke nicht definiert.
4.3.1	Visualisation of the measurement data.....	Fehler! Textmarke nicht definiert.
4.3.2	Definition of the steps	Fehler! Textmarke nicht definiert.
4.3.3	Printing	Fehler! Textmarke nicht definiert.
4.4	The Report (Step 4)	Fehler! Textmarke nicht definiert.
4.4.1	Explanation of the report.....	Fehler! Textmarke nicht definiert.
4.4.2	Data rating	Fehler! Textmarke nicht definiert.
5	FEATURES OF REPORT.....	21
6	APPENDIX	22
6.1	Minimum computer equipment.....	22
6.2	Software-Installation of BalanceLab.....	22
6.3	Hardware Installation	25
6.4	Backup of all measurement data	26
6.5	Import and export of data.....	27
6.5.1	Import and export of individual measurements	27
6.5.2	Export ASCII and APD.....	28
6.5.2.1	Export ASCII.....	30
6.5.2.2	Export APD.....	32
6.6	Notes for measuring principle	34
6.7	Safety instructions.....	35

1 Introduction

The BalanceLab program for analyzing the distribution of static and dynamic forces runs under Windows XP and Windows Vista.

Please find the instructions on the computer hardware requirements in the appendix.

1.1 Layout of the program



The program includes a data base (1) in which projects, players, and individual measurements are filed. On the basis of these data, the measurements (2) are run. Afterwards the viewer (3) offers the possibility to review the recorded data again as well as in slow-motion or time lapse. Different calculated data can easily be accessed with a few clicks. This is done in the report function (4). Moreover, the program includes a measuring data export option.

2 Preparing the database (step 1)



Note

Please first install the BalanceLab-Software before you plug in the USB cable of the platform.
Information about the installation is given in the appendix chapter 6.2

After the installation of the BalanceLab-platform please connect the USB cable between the PDMS platform and the computer. Double-click the BalanceLab icon or launch the program with "Start" - "Programme" - "SAM" - "BalanceLab". The BalanceLab program will now open the database screen and you can prepare your measurement.

2.1 Instructions

The screenshot shows the 'Database' application window. It features a 'Project' dropdown menu set to 'Beispieldaten', with 'New', 'Rename', and 'Delete' buttons. Below this is a 'Player' list containing 'Beispiel' and a 'Records' list with four entries: '02-03-2006 Stance Test', '02-03-2006 Dynamic Test rechter Fuß (4x)', '02-03-2006 Dynamic Test rechts und links', and '02-03-2006 Stance Test 2'. A right-hand sidebar contains 'Setup', 'Measure', 'View', 'Report', and 'Exit' buttons. A bottom toolbar includes 'New', 'Edit', 'Delete', 'Copy...', 'Rename', 'Delete', and 'Copy...' buttons. Three numbered callout boxes provide instructions: '1. Create a project' points to the 'New' button; '2. Create a new entry with player's data' points to the 'New' button in the bottom toolbar; '3. Click "Measure"' points to the 'Measure' button in the sidebar.

2.2 General information

The BalanceLab database system has three levels for measurement data management:

The highest level is the project level. Here various projects, patient groups or program users can be assigned. Players` names are entered in the next level ("Players"). Measurement data files are listed in the last level ("Records").

The measurement file names are automatically assigned by the program according to the date of examination and a consecutive number (e.g. 02-03-2006 Stance Test 2). It is possible to rename data files.

By clicking on "New", "Delete", "Rename" or "Edit" projects, players or records may be added, deleted or renamed.

2.3 Create new players


By selecting the functions "Player" "New" the following dialog appears for entering the player`s data:


When entering a player's name, the sex and the birth date as well as a player specific code can be entered. The birth date is entered under "Born" in the following format: DD.MM.YYYY (the separation can be done with . / or space, e.g. 21/03/06 and the program will automatically replace them with dots: 21.03.2006). The field "Comments" can be used for findings, memos, etc.


3 Analyzing the pressure distribution

The BalanceLab system analyses the static pressure distribution in real-time either two- or three-dimensional. Furthermore it is possible to record an unlimited number of measurements one after the other and to compare two different records in one report.

3.1 Instruction for the pressure analysis

To start the measurement please select "Stance Test" and click on  .


The measurement takes about 20 seconds. To adjust the measurement duration, please click  under "Configuration". Here the presetting of 20 seconds can be changed.

	<p>Please do not stand on the platform, when the Start-Button is pressed. Please wait until the hardware initialization has finished and consider the given instructions.</p>
<p>NOTE</p>	



First the platform has to be calibrated. For this, please leave the platform and press "Enter".


Afterwards following screen will be shown



Window A displays the pressure distribution under the feet 2D or 3D. The color scale on the left representing the force in N/cm² acting on the individual plate sensors. The pink line is the connecting line of the three CoPs (Center of pressure). The middle cross is the center of pressure of the contact area. Left and right points are the CoPs of the left and right contact area.


Window B shows the force curves in relation to the measuring time. The vertical reaction force of the left and right forefoot is shown at the top and underneath the reaction force of the left and right backfoot.

In *Window C* the percentage values and the red and green bars indicate the force-distribution between the left and right foot and the distribution between left and right front and back foot .




Window D shows the various instructions for the respective measuring process and a scale which indicates the time of measurement.


3.2 Execution of measurement (step 2)



Firstly the player stands on the platform. For the correct positioning please consider the label on the platform.




The red lines separate the left and right body side as well as the fore- and backfoot.


For correction, please press the  button. Afterwards the lines can be shifted. This step will be explained more detailed in chapter 3.3.2.

Please press the red record button  to start the measurement. Alternatively the “Enter” key can be used.

The measurement is stopped automatically after the selected measurement duration.

To abort the measurement at an earlier time, use the Stop button  or press the “Enter” key.

When the record is finished following dialog box appears.





Here the name of the record can be input. After confirming with “OK” the measurement will be saved and a new measurement can be performed.


The number of recorded measurements is not limited.

By clicking on “Cancel” the last performed measurement will not be saved and the main database is shown.

3.3 Preparation for the report in the „View” (step 3)


The “View” program section is achieved via the database by clicking on the  button. Following screen appears:



The View screen shows first or after activating the button  the average load distribution under the feet.






The pink line is the connecting line of the three CoPs (Center of pressure). The middle cross is the center of pressure of the whole body. Left and right points are the CoPs of the left and right body side.


The distance which the CoPs cover during the measurement are drawn as white curves.

With the print button  the current screen view can be print out. In doing so it should be considered that the print job is send to the default printer of the computer.


3.3.1 Visualization of the measurement data




Several features are available for the visualisation of the measurement data:

- Automatic play-back of the measurement by clicking the Play button  or pressing "Enter". The measurement is played-back and repeated until the Pause button  (or the "Enter" key) is pressed.
- Press the  button to see the current force distribution. By shifting the dashed line (with the mouse or the right and left arrow keys) in the force/time diagram the measurement can be viewed at different points of time.
- With the "zoom" buttons   you can zoom in/zoom out the presentation of the platform to the left or the force/time curve to the right.
- By clicking with the left mouse button, holding the shift key and clicking the left mouse button again, an interval in the force/time diagramm can be selected (blue background). Now play-back or zooming of the interval is possible.
- The scaling of the force/time diagram can be modified. For doing so the mouse must be moved on the force axis of the force/time diagram so that the following symbol appears:




By keep holding the left mouse button and moving left and right the scaling can be adapted. In the same way the position of the base line can be changed by moving up or down.

- The view of the platform can be changed in any way. First you have to activate the button . For changing the projection following features can be used:

	<p>Enables the scrolling of perspective views with the mouse movement (hold the left mouse button).</p>
	<p>Enables the zooming of views by moving the mouse (hold the left mouse button).</p>
	<p>Enables the rotation of perspective views. To achieve this, move the mouse (hold the left mouse button).</p>

3.3.2 Preparation

For preparing the report the software needs to know, where the left and right foot and the forefoot and backfoot is positioned. This can be defined by clicking the  button (the lines turn yellow) and move the vertical and horizontal lines that way, that the necessary information is correct.


To get to the report, please leave the viewer by clicking the Exit-Button. Please confirm the question for “Save Changes” with Yes.


3.4 The report (Step 3)



To get to the report please click the button, which is available in the data base.

3.4.1 Explanation of the report






3.4.2 Comparison of two records


In order to compare two different measurements, please proceed as follows:

- Mark the first record
- Mark the second record of equal type (stance test or dynamic test). To mark the second record keep the 'Ctrl' key pressed and click on the second record.
- Click the 'Report' button in the database



The report shows now colored and grey diagrams. The colored diagrams represent the first and the grey diagrams represent the second measurement.

3.4.3 Data Rating



The image above indicates the following:

The load distribution between left and right foot is expressed as a percentage.

A perfect distribution is an equal distribution of 50 % for the right and left contact area.

Ideally, the connecting line of left and right CoP (pink line) runs parallel to the connection of the heel contacts.


By the load distribution between backfoot and forefoot the forefoot should be loaded with approx. 1/3 (33%) and the backfoot with approx. 2/3 (66%).

The maximum load should not exceed 15 N/cm².

For sensitive feet (e.g. diabetics, etc.) the pressure should be less than 11 N/cm² in order to avoid any damage to the soles of the feet.

These values apply to a “normal” load distribution during upright standing without shoes.

4 Features of Report

To change the layout of the report click on the “Edit Layout“ button .

- **Size and position of frames**

Size and position of the graphic presentations can be modified by clicking on the frames and dragging with the mouse.

- **Header and footer**

Under “Edit” – “Header” / “Footer” the top or bottom line can be changed.

- **Saving**

New settings are stored with "Save Template" under "File".

- **Zoom buttons** 

With the zoom buttons, the actual size of the report can be maximized or minimized.


- **Viewing the whole page** 

By clicking the button for viewing the whole page, the report can be shown on the screen as it will be printed. With the zoom buttons, the actual view can be maximized or minimized. You can view "Two Pages" at the same time or change to the "Next Page" or "Prev Page" by clicking the buttons.

With "Close" the side view is finished.

With “Print...” the report can be printed. Please note that for this operation, the driver of the printer must be installed.

- **Exit** 

Leave the Report screen by clicking on the “Exit” button or on  in the upper right corner of the window.

5 Appendix

5.1 Minimum computer equipment

- Intel Pentium IV or better
- 1024 MB RAM
- VGA graphics card with resolution 1024 x 768 and 24Bit True color
- Color monitor with minimum resolution 1024 x 768
- Windows XP or Windows Vista
- 200 MB free disk space on drive C:
- min. one free USB port

5.2 Software-Installation of BalanceLab



NOTE


Attention: Please do not connect the USB cable between the PDMS platform and the computer before the software has been installed on your computer!

To install BalanceLab start Windows and place the CD in the computer. The installation runs automatically.

Following window opens:



By clicking on the "Next" button the installation will start.



Accept the license agreement and click "Next" to continue.



By default the program will be installed to the folder " SAM BalanceLab 1.0 into your Windows programs folder. You can change the installation location if you like.

In the case of installing the software as an update, the existing measuring data will be included in the update automatically

After successful installation an icon with the name "SAM BalanceLab" automatically appears on the desktop.



Now at this stage (but not beforehand), please connect the platform to the computer using the enclosed USB cable !


NOTE

After this, the USB driver will be installed automatically. Please ignore any error messages of the Windows operating system that might appear during the installation. Be sure to press "Continue" at any warning message, never press "Cancel".

5.3 Hardware Installation


Normally you would not need to change anything here. The default values are set for operation with USB connection.

Select „Setup“  **Setup** in the database. You will reach the hardware setup:



The interface mode (serial or USB) can be changed with „Others“.

To select the measurement device (USB interfaces) or the port (e.g. COM1) use „Settings“.



For information about your measurement device please consider the table at the platform side next to the cable connections.

NOTE

This table shows which settings have to be confirmed before measuring:


connection cable	Other...	Settings...
USB	SDK Driver	FDM
Serial	Serial	choosing the port
	SDK Driver	„PDM device on RS232“ and choosing the port with „Setup“

5.4 Backup of all measurement data

In order to avoid loss of data due to a technical defect of the hard disk a backup of measurement data in regular intervals is recommended.

The program path is set to C:\Programme\SAM\BalanceLab on standard installation.


In this directory, the whole data are arranged in the sub-directory „Data“.



In order to back up all the measurement data, this sub-directory is copied on to a second hard disk, tape drives, DVD or CD. Special backup software is available for the tape drives, which is normally supplied with the drive. Backup on a second hard disk can be performed with the help of the Windows Explorer.

5.5 Import and export of data

The database has a „Copy“ function to import or export players with their measurement data or just a specific measurement. Furthermore it is possible to export all data in different formats.




For single records following features are available:


- „External Records“
- „APD Export“
- „ASCII Export“

5.5.1 Import and export of individual measurements

By clicking the „Copy“ button from “Players” or by clicking the “Copy” button from “Records” and selecting “External Records” the following dialog box appears:



After confirmation the next dialog box appears:



The left side represents the existing, local database. On the right-hand side any database can be selected which is on a floppy disk, hard disk or on other data carrier. The background of the respective selected record is blue. By pressing the „<<“ or „>>“ button the selected data are imported in the local database or are exported in the chosen directory.

Whereas the double-arrow to the left << is for the import and the double-arrow to the right >> is for the export of the measurement data.


5.5.2 Export ASCII and APD

Similar to the normal export (chapter 6.5.1) a target folder for saving the exported files has to be selected first.


In the following chapter the structure of the ASCII- and APD-files are explained.

For better understanding of both export results, following pictures show the position of the used coordinate system in relation to the PDMS platform.

View static measurement:



View dynamic measurement :



5.5.2.1 Export ASCII

By exporting ASCII data a text file is created. The structure of the ASCII-File is built up as follows (e.g. dynamic analysis):

```

1 Person Renate Beispielmann
  Sex F
  Born 26/03/1983
  Code
  LegLength 84
2 Record 06-03-2006 Dynamic Test mit Schuh
  Application DynPlat
  Creation date 06/03/2006 13:52:37
3 Exercises {
  # Name Start,sec Length,sec Start time
  1 Right 0.00 0.90
  2 Left 0.90 0.88
  3 Right 1.78 0.85
  4 Left 2.63 0.88
  }
4 Frequency,Hz 60.000
  Count 211
5 Time,ms Force,N Max pressure,N/cm^2
  Calibr
  0 40.0 12.5
  17 59.0 18.0
  33 98.5 24.0
  50 115.5 30.0
  67 149.5 36.0
  83 183.5 35.5
  100 234.0 37.5
  117 287.0 39.5
  133 352.5 29.0
  150 457.0 19.0
  167 528.0 24.5
  183 574.0 30.0
  200 578.5 30.0
  
```


The main part of the file consists of frames of each recorded sample. In the frame the force value of each sensor is shown. The example shows, that 211 samples were recorded. The number of samples depends on the frequency and the duration of the measurement (here 60 Hz and 3,51sec)

7	Title of the chart, consecutive numbering of the charts, position of the center of pressure CoP
8	Matrix of recorded forces for each sensor. Columns= x coordinate, lines =y-coordinate
9	Continuation of the chart

5.5.2.2 Export APD

In case of an export to APD format is performed, files with the extension „.apd“ are created. If a dynamic measurement is exported to APD-format and more than one roll-off analysis was performed, each roll-off analysis is exported to one single APD file.


```

[General]
Description=pedcad Importfile
Source=<<zebris>>
Date=06/03/2006
Time=13.52.37
UnitDistance=mm
UnitPressure=N/cm2

[Customer]
FirstName=Renate
LastName=Beispielmann
Sex=F
FootSide=Right

[Technical]
MaxSensorsX=47
MaxSensorsY=32
StartSensX=15
StartSensY=10
SensCountX=29
SensCountY=11
LDistX=10
LDistY=10

```



1	general information about the measurement, file
2	player data (name, sex, foot side)
3	technical platform and sensor data: <ul style="list-style-type: none"> • maximum number of sensors in x and y direction • the lowest sensor number, seen from the origin, which was loaded • number of heavily loaded sensors in x and y direction • sensor size

4 Chart of recorded forces for each sensor (only the part, which was loaded). Columns= x coordinate, lines =y-coordinate. The value „-1“ stands for an unloaded sensor. Please consider the position of the coordinate system.

5.6 Notes for measuring principle



NOTE

The measuring system must be checked at regular intervals to ensure that the measuring function is operating correctly.

Notes

If invasive measures are to be taken, the measuring system must only be used as additional expert evidence. On no account can or may invasive interventions or measures potentially dangerous to the player be undertaken solely on the basis of the measuring results.

In the event of recognizable damage to the sensor area of the measuring platform, no further measurements must be undertaken.

After the performance of a zero measurement, no measured values must be displayed in the unloaded state. Furthermore, the pressure distribution images must be examined regularly for atypical measuring patterns. These mainly concern column or row shaped measuring patterns distinguishable from the surrounding values. If any of the faults described occur or in the event of doubt, the manufacturer must be contacted.

Calibration

The measuring accuracy of the sensors for the pressure distribution measurement must be checked from time to time using a specific force.

This can be done by a user of known body weight standing on the platform with one foot. Taking into account the acceleration due to gravity, the platform must indicate any not fully loaded marginal sensors and the measuring tolerance of the approximate body weight. A false indication means that re-calibration must be carried out by the manufacturer. In any event, the measuring platform must be checked and re-calibrated every two years by the manufacturer.

5.7 Safety instructions



CAUTION

General safety instructions for operating the system

- Never use the device in a damp place or where it could possibly be sprayed with fluids. Any fluid penetrating the device can cause a fire, electric shock or other serious accidents.
- Never pour any fluid over the system or its components.
- Do not place the system on an unstable surface.
- Do not install the system above or near any heating.
- Only operate the system using the stipulated mains voltage and the supply unit that has been approved for medical products and is included in the delivery by the manufacturer.
- Lay the mains cable such that no-one can trip over it and check regularly that it is not damaged.
- Never insert any objects in the system components.
- Do not attempt at servicing the system in any manner other than that described in this manual. If the cover is removed it is possible that you could be exposed to lethal voltages or other hazards.
- Should one of the situations in the following list occur, please do not fail to contact the manufacturer:
 - a) The mains cable or mains plug has been damaged
 - b) Fluid has been spilled over the main device or other system components
 - c) The system fails to function properly although the operating instructions have been adhered to
 - d) The platform or other system components has/have been dropped or the housing has been damaged



CAUTION

Safety instructions for applying the PDM-S Platform on the player results.

- Before using the device for the first time, please read the operating manual and safety instructions carefully. If there is any uncertainty about applying the device correctly, or about its correct mode of function, the Instruction Manual must always be consulted. Therefore the Instruction Manual should always be kept within easy reach of the product.

- The application and operation of the system and also the evaluation of the measuring data and their interpretation may only be carried out by trained qualified personnel. The manufacturer assumes no liability for any injury to persons, damage to property, or loss of data due to improper use of the software, the device or its component parts.
- The players' data and measuring data may only be copied, moved, or deleted using the database function provided by the SAM application programs. In the case of data being changed intentionally without using the database functions, the user alone bears the full risk.
- Before beginning each measurement, check to ensure that the correct choice and correct position of the transmitters or application aids have been carried out.
- An injury hazard exists for the player from the cables. Please observe here the special instructions in the application software manuals and take particular care not to allow any children or mentally retarded persons to go near the device without supervision.
- Should there be any damage to the device or component parts, they should be returned to the manufacturer for a safety check. It is forbidden to continue using them, as severe damage and serious injuries - even lethal injuries - may result. The manufacturer must always be contacted in all cases of fault or doubt.
- We also point out that by making changes to this certified device or its accessories your legal right to operate it will be nullified.
- Should any invasive measures be taken, the measuring system may only be applied as a supplementary means for expert evaluation. On no account must or may invasive measures or measures endangering the player be carried out solely on the basis of the measuring results.
- The platform must be set up on a non-slip base, or built into a catwalk, in order to rule out any danger to the player due to the platform slipping.
- The servicing, repair and re-adjustment must only be carried out by authorized, qualified specialist. Only store and transport in the original packing supplied by the manufacturer.